OSTEOLIGN

The Premier Drill Guide for Arthroscopic Transosseous Rotator Cuff Tear Repair

The Problem and Need

- Suture Anchors = \$\$\$\$ \diamond
- < 90° critical angle=
- **REVISION SURGERY!**

CLEAT Engineering's Solution:

- Bypass the critical angle ENTIRELY \diamond
- **Deliver STRAIGHT, uniform 60° bone** \diamond tunnel
- 60° has shown to deliver tunnels \diamond with adequate bone roof while bypassing both the axillary nerve



Why OSTEOLIGN?

Straight bone tunnel deliverv

Mitigates Bone-Suture Cut out

Prevents Revision Surgery

- **ELIMINATES NEED FOR SUTURE** \diamond **ANCHORS**
- Hospitals save \$956/surgery on * average
- \$900 estimated ROI/surgery *
- Reusable: 50+ Surgeries
- **Compatible with currently** available arthroscopic tools

The Market:

- The soft tissue repair market is expected to increase by 6% from 2017-2025 reaching \$9.7 billion in the U.S.
- 250,000 RTCT repairs performed annually *
- 78% of all RTCT repairs are arthroscopic *
- \$3.3 Billion annual burden produced due to * **RTCT** repair











Contact CLEAT Engineering: Colleen Martin 843-697-3619 | collee5@g.clemson.edu Liz Donley 412-478-5147 | edonley@g.clemson.edu Erin Lester 774-291-0407 | emleste@g.clemson.edu Andrew Losier 551-427-1326 | alosier@g.clemson.edu Taylor Sanders 630-460-1226 | tpsande@g.clemson.edu



